

## L Laboratory Furnaces Technical Specifications

Furnace type	Tmax °C	Recommended operating temperature in °C	Volume in L	External dimensions (w × h × d) mm	Internal dimensions (w × h × d) mm	Input in kW	Weight in kg	Protection*	Voltage V	Max load capacity of bottom in kg
L 03/12	1200	700 - 1150	3	380×440×400	180×100×140	1.2	20	16/1	230	4
L 05/12	1200	700 - 1150	5	430×470×430	230×130×170	2.4	26	16/1	230	6
L 09/12	1200	700 - 1150	9	430×505×500	230×170×240	3.0	30	16/1	230	6
L 15/12	1200	700 - 1150	15	450×505×600	250×170×340	3,5	40	16/1	230	6

\*Protection can vary depending on accessories added (at additional cost).

Subject to technical changes.