

LE Laboratory Furnaces Technical Specifications

Furnace type	Tmax °C	Recommended operating temperature in °C	Volume in L	External dimensions (w × h × d) mm	Internal dimensions (w × h × d) mm	Input in kW	Weight in kg	Protection*	Voltage V	Max load capacity of bottom in kg
LE 05/11	1100	700 - 1050	5	470×330×455	170×130×230	1,8	20	16/1	230	6
LE 09/11	1100	700 - 1050	9,4	490×370×515	190×170×290	2,3	24	16/1	230	6
LE 15/11	1100	700 - 1050	14,4	550×370×565	250×170×340	3,0	28	16/1	230	6

*Protection can vary depending on accessories added (at additional cost).

Subject to technical changes.